

Science Safety, Inc. 7345 164th Ave NE STE i145-1283 Redmond, WA 98052 https://sciencesafety.com

Incompatible Chemical Typical Table

Chemical	Incompatible With These Chemicals
Acetic acid	Chromic acid, nitric acid, hydroxyl compounds, ethylene glycol, perchloric acid, peroxides, permanganates
Acetic anhydride	Hydroxyl-containing compounds such as ethylene glycol, perchloric acid
Acetone	Concentrated nitric and sulfuric acid mixtures, hydrogen peroxide
Acetylene	Chlorine, bromine, copper, fluorine, silver, mercury
Alkali and alkaline earth metals such as powdered magnesium, sodium, potassium	Water, carbon tetrachloride or other chlorinated hydrocarbons, carbon dioxide, halogens
Ammonia (anhydrous)	Mercury, halogens, calcium hypochlorite, hydrofluoric acid
Ammonium nitrate	Acids, metal powders, flammable liquids, chlorates, nitrites, sulfur, finely divided organic or combustible materials
Aniline	Nitric acid, hydrogen peroxide
Arsenical materials	Any oxidizer agent
Azides	Acids, heavy metals and their salts, oxidizing agents
Calcium oxide	Water
Carbon, activated	All oxidizing agents, calcium hypochlorite
Carbon tetrachloride	Sodium
Chlorates	Ammonium salts, acids, metal powders, sulfur, finely divided organic or combustible material
Chlorine dioxide	Ammonia, methane, phosphine, hydrogen sulfide
Chromic acid and chromium trioxide	Acetic acid, alcohol, camphor, glycerol, naphthalene, flammable liquids in general
Copper	Acetylene, hydrogen peroxide
Cumene hydroperoxide	Acids (organic or inorganic)
Cyanides	Acids
Flammable liquids	Ammonium nitrate, chromic acid, hydrogen peroxide, nitric acid, sodium peroxide, halogens, other oxidizing agents
Fluorine	All other chemicals
Hydrides	Water
Hydrocarbons (e.g., butane, propane, benzene)	Fluorine, chlorine, bromine, chromic acid, peroxides

Hydrocyanic acid	Nitric acid, alkalis
Hydrofluoric acid (anhydrous)	Ammonia (aqueous or anhydrous)
Hydrogen peroxide	Copper, chromium, iron, most metals or their salts, any flammable liquid (i.e., alcohols, acetone), combustible materials, aniline, nitromethane
Hydrogen sulfide	Fuming nitric acid, oxidizing gases
Hypochlorites	Acids, activated carbon
lodine	Acetylene, ammonia (aqueous or anhydrous), hydrogen
Mercury	Acetylene, fulminic acid, ammonia
Metal hydrides	Acids, water
Nitrates	Acids
Nitric acid (concentrated)	Acetic acid, acetone, alcohol, aniline, chromic acid, hydrocyanic acid, hydrogen sulfide, flammable liquids, flammable gases, copper, brass, any heavy metals
Nitrites	Acids
Nitroparaffins	Inorganic bases, amines
Oxalic acid	Mercury and silver and their salts
Oxygen	Oils, grease, hydrogen; flammable liquids, solids, or gases
Perchloric acid	Acetic anhydride, alcohol, bismuth, paper, wood, grease, oils
Permanganates	Concentrated sulfuric acid, glycerol, ethylene glycol, benzaldehyde
Peroxides, organic	Acids (organic or mineral), avoid friction, store cold
Phosphorus, white	Air, oxygen, alkalis, reducing agents
Potassium	Carbon tetrachloride, carbon dioxide, water
Potassium chlorate	Sulfuric and other acids, ammonium salts, metal powders, sulfur, finely divided organics, combustibles
Potassium perchlorate (see also chlorates)	Sulfuric and other acids
Potassium permanganate	Glycerol, ethylene glycol, benzaldehyde, sulfuric acid
Silver and silver salts	Acetylene, oxalic acid, tartaric acid, ammonium compounds, fulminic acid
Sodium	Carbon tetrachloride, carbon dioxide, other chlorinated hydrocarbons, water
Sodium nitrate	Ammonium nitrate and other ammonium salts
Sodium peroxide	Ethyl or methyl alcohol, glacial acetic acid, acetic anhydride, benzaldehyde, carbon disulfide glycerin, ethylene glycol, ethyl acetate, methyl acetate, furfural
Sulfides	Acids
Sulfuric acid	Chlorates, perchlorates, permanganates

Summary has been adapted from <i>Prudent Practices in the Laboratory: Handling and Disposal of Chemicals</i> , National Research Council, 1995